



THE CRYO LC



CRYO LC

INTRODUCTION

The CRYO LC is a technologically advanced localized cryotherapy device designed for **effective, convenient, and safe targeted cryotherapy** applications.

This device is equipped with advanced features for safe and ergonomic everyday use, such as **laser distance measurement** and **skin temperature reading**. Its compact size and ergonomic design make it an ideal device for everyday use.

The device uses liquid nitrogen to achieve cryogenic temperatures offering four operating modes ranging from **-80°C to -160°C / -110°F to -250°F**. As a result of the flexible temperature settings, a wide range of applications is possible.

The CRYO LC is an excellent addition to **physiotherapy clinics, gyms, sports clubs, spas, and beauty salons**.



INTRODUCTION

LOCALIZED CRYOTHERAPY

Localized cryotherapy is the application of cryogenic temperatures onto a specific area of the body. This technique allows for the **targeted treatment of specific parts of the body**, such as the neck, elbow, knee, ankle, face, and others.

The session **only takes minutes** to complete and leaves clients with benefits that far exceed what traditional icing can provide.

What are the applications for localized cryotherapy?

Localized cryotherapy is commonly used for pain management, injury recovery, and aesthetic purposes. It is used most frequently in **physiotherapy and aesthetics**.



**Muscle
Recovery**



**Health
& Wellness**



**Pain
Management**



**Beauty
& Anti-aging**



**Weight
Loss**

REPORTED BENEFITS

PHYSIOTHERAPY

- Pain relief
- Reduction of swelling
- Reduction of inflammation
- Improved blood circulation
- Faster recovery

AESTHETICS

- Reduction of cellulite
- Increase in collagen production
- Improvement of skin texture
- Enhancement of youthful appearance
- Faster recovery following cosmetic procedures

CRYO LC

HIGHLIGHTS

Operator Touchscreen

Simple & intuitive session management

50 L Nitrogen vessel

Up to 40 sessions capacity

Custom LED lighting

12 ambient color options

4 Session modes

Custom intensity settings

Skin temperature display

Real-time skin temperature reading

Laser distance control

Laser guidance for optimal distance

Tank level indication

Real-time nitrogen level indication



CRYO LC

FEATURES

FEATURE	DETAIL	INCLUDED
OPERATION	Compact dimensions & light weight Ergonomic, light-weight handle Large, user-friendly touchscreen control panel	✓
SAFETY FEATURES	Laser distance indication Skin temperature reading Audible notifications	✓
CUSTOMIZATION	2 design vessel covers – titanium & white 12 LED Light color options	✓
NITROGEN VESSEL	50L vessel (approx. 400 minutes of operation)	✓
SIMPLE REFILLING	Refilling is possible without dismounting the vessel	✓
NO ASSEMBLY	Comes fully assembled	✓

DIMENSIONS



Weight: 30 kg / 66 lb (empty)
70 kg / 154 lb (full)



RECOVERY FOR EVERYONE

Contact us

info@cryoniq.com

Visit us

www.cryoniq.com

